

Annex M - Equipment List - Springhill Institution / Annexe M - Liste d'équipements - Établissement Springhill

	A	B	C	D	E	F	G
1	Item	~~Equipment Code~~Code d'équipement~~	~~Institution~~Établissement~~	~~Location~~Localisation~~	~~Manufacturer~~Nom du fabricant~~	~~Model~~Modèle~~	~~Serial Number~~Numéro de série~~
2	AccuVein HF580	TBD	Springhill Institution	TBD			
3	AccuVein HF580	TBD	Springhill Institution	TBD			
4	Air Filter for Air Purifiers		Springhill Institution	Dentist	Ausin Air Systems Canada Inc		
5	Air Purifier HM400		Springhill Institution	Dentist	Austin Air Systems Canada Inc	HM400	601425181
6	AMALGAMATOR-DENTAL	BMD 10538	Springhill Institution	Dentist	Ivoclar Vivadent Inc	15477 Silamat 56	108955
7	AMALGAMATOR-DENTAL	BMD 11644	Springhill Institution	dentist	Ultra Mat 2		71903
8	ASPIRATOR	BMD 7869	Springhill Institution	Treatment	Devil Bless	7305	PD235889
9	ASPIRATOR-DENTAL	BMD 5351	Springhill Institution	Treatment	Gomco	4005	20071218026
10	CAMERA-CONTROLLER-LIGHT-SOURCE	BMD 11646	Springhill Institution	Back Office	SomeTech Inc	DCSM(XGA)	DCSN-11B-002
11	CAMERA-VIDEO	BMD 11647	Springhill Institution	Back office	Sometech Inc	DCSN	DCSN-11C-117
12	CAMERA-VIDEO	BMD 11648	Springhill Institution	Back Office	SomeTech	DCS-104(NTSC)	DCSN-11B-058
13	CART-DENTAL	BMD 7855	Springhill Institution	Dentist	ADEC		
14	CHAIR-DENTAL	BMD 7846	Springhill Institution	dentist	Belmont	X-Calibur	73203603
15	CHAIR-OPHTHALMIC	BMD 5349	Springhill Institution	Optometrist	Topcon	238717	222466
16	COVID-19 Swab Testing Machine		Springhill Institution	Health Services	Abbott	Abbott ID Now	4502D91C
17	COVID-19 Swab Testing Machine		Springhill Institution	RRC	Abbott	Abbott ID Now	7362E61C
18	Dental - Drill handpiece		Springhill Institution	Dentist	Kavo	LUX E675 L	
19	Dental Drill Handpiece		Springhill Institution	Dentist	Kavo	LUX E675L	
20	DENTAL-UNIT-ULTRASONIC	BMD 16065	Springhill Institution	Dentist	L&R Medical	Quantrex 140	3021822
21	DOPPLER-FETAL-HR	BMD 15960	Springhill Institution	NP Office	Summit Doppler System	L250	L250
22	ELECTROCARDIOGRAPH	BMD 16025	Springhill Institution	Doctor Office	GE	MAC5500 HD	SKJ20090008PA
23	ELECTROCARDIOGRAPH	BMD 7856	Springhill Institution	Treatment	GE	MAC5500	SCD10087696PA
24	HYDROCOLLATOR-HOT-PACK	BMD 5352	Springhill Institution	Physio	Chattanooga Canada	Hydrocollator -2102C	2328C
25	IMAGER-COMPUTER-RADIOGRAPHY	BMD 10435	Springhill Institution	X-Ray	AGFA	CF30X	1380
26	INTERFERENTIAL-ULTRASOUND	BMD 5333	Springhill Institution	Physio	Chattanooga Canada	2762CC Intellect Advanced	4295
27	KERATOMETER	BMD 11641	Springhill Institution	Optometrist	Topcon	238714	01074G
28	LAMP-SLIT	BMD 11643	Springhill Institution	Optometrist	Topcon	238725	203100017B
29	LENSOMETER	BMD 7543	Springhill Institution		Bausch & Lomb	712655	512045
30	LENSOMETER		Springhill Institution	Health Services	Right Med	Model 1 S4OPTIK	
31	LIGHT-EXAMINING	BMD 15906	Springhill Institution	Treatment	Welch Allyn	2023	N/A
32	LIGHT-EXAMINING	BMD 15910	Springhill Institution	NP office	Welch Allyn	GS300	
33	LIGHT-EXAMINING	BMD 7863	Springhill Institution	Treatment	Welch Allyn	44300/LS135	
34	LIGHT-EXAMINING	BMD 7863	Springhill Institution	Doctor Office	Welch Allyn	3195 - LS135	N/A
35	LIGHT-EXAMINING (CEILING)	BMD 7852	Springhill Institution	dentist	Belmont 501T	105000 Dental Light	H71903
36	LIGHT-HEAD	BMD 11640	Springhill Institution	Optometrist	Heine	268299	0/11848
37	LIGHT-ULTRAVIOLET	BMD 15200	Springhill Institution	Dentist	Kerr	Demi Plus	
38	LIGHT-ULTRAVIOLET	BMD 7853	Springhill Institution	dentist	Bobcat	506 QHL75	506-30012
39	MONITOR-NIBP	BMD 15908	Springhill Institution	Treatment	LifeSource	UA-767FAM	7190700142
40	MONITOR-NIBP-SPO2-TEMP	BDM-100159	Springhill Institution	Health Services	Welch Allyn	Spot Vital Signs 4400	2110009420421
41	MONITOR-NIBP-SPO2-TEMP	BMD 10430	Springhill Institution		Welch Allan	42NTB	
42	MONITOR-NIBP-SPO2-TEMP	BMD 7978	Springhill Institution	Doctor Office	Welch Allyn	42MTB Spot Vital	201103450
43	MONITOR-NIBP-SPO2-TEMP	BMD 8084	Springhill Institution	RRC - Nursing Office	Welch Allyn	420 Series	201216120
44	MONITOR-NIBP-SPO2-TEMP		Springhill Institution	Telemedicine room	Welch Allyn	42NTB	201739756
45	MONITOR-NIBP-SPO2-TEMP	BMD 8083	Springhill Institution	Doctor Office	Welch Allyn	420 Series	201216780

Annex M - Equipment List - Springhill Institution / Annexe M - Liste d'équipements - Établissement Springhill

	A	B	C	D	E	F	G
1	Item	Equipment Code Code d'équipement	Institution Établissement	Location Localisation	Manufacturer Nom du fabricant	Model Modèle	Serial Number Numéro de série
46	MONITOR-NIBP-SPO2-TEMP	BMD 10518	Springhill Institution	Treatment	Welch Allyn	Vital Signs LXI	20150100638
47	MONITOR-NIBP-SPO2-TEMP	BMD 10528	Springhill Institution	Treatment Room/Portable	Welch Allyn	42NTB	201403441
48	MONITOR-NIBP-SPO2-TEMP	BMD 10529	Springhill Institution	A&D	Welch Allyn	0297	201403328
49	MONITOR-NIBP-SPO2-TEMP		Springhill Institution	Health Services	Welch Allyn	WA-42NtB-E1	2019-08614
50	MONITOR-NIBP-SPO2-TEMP	BMD 16300	Springhill Institution	NP Office	Welch Allyn	Connex VSM 6400	
51	MONITOR-SPO2	BMD 10526	Springhill Institution	RRC	GE	Ohmeda TuffSat	FCB08110053SA
52	MONITOR-SPO2	BMD-15909	Springhill Institution	Treatment	G E Healthcare	Ohmeda Tuffstat	
53	MONITOR-SPO2	BMD 15907	Springhill Institution	Treatment	GE	Ohmeda Tuffstat	FCB08110036SA
54	MONITOR-SPO2	BMD 5337	Springhill Institution	Treatment	Nellcor Puritan Bennett	NPB40	G00842991
55	MONITOR-SPO2	BMD 7857	Springhill Institution			365367	G00842987
56	MONITOR-SPO2	BMD 8080	Springhill Institution	Doctor Office	Nellcor	N-65P	P12827521
57	MONITOR-TEMP	BDM 10523	Springhill Institution	Treatment	Welch Allyn	Braun Pro4000	23715K01237
58	MONITOR-TEMP	BDM 10525	Springhill Institution	Treatment	Welch Allyn	Braun Pro 4000	32414K03641
59	MONITOR-TEMP	BMD 5345	Springhill Institution		Welch Allyn	679 Suretemp	3137764
60	MONITOR-TEMP TYMPANIC	BMD 10520	Springhill Institution	Treatment	BRAUN	6022	11111140252
61	MONITOR-TEMP-TYMPANIC	BMD 10425	Springhill Institution		Welch Allan	6022	
62	MONITOR-TEMP-TYMPANIC	BMD 10426	Springhill Institution		Welch Allan	6021	
63	OPHTHALMIC-INSTRUMENT-STAND	BMD 5350	Springhill Institution	Optometrist	Topcon	238721	350236
64	OPHTHALMOSCOPE	BMD 11645	Springhill Institution	Doctor Office	Welch Allyn Canada	767	
65	OPHTHALMOSCOPE	BMD 5346	Springhill Institution		Welch Allyn	767	
66	OPHTHALMOSCOPE	BMD 7871	Springhill Institution	Doctor Office	Welch Allyn	767	
67	OPHTHALMOSCOPE (Integrated Diagnostic System)	BMD 15905	Springhill Institution	RRC	Welch Allyn	GS 777	
68	OPHTHALMOSCOPE (Integrated Diagnostic System)	BMD 16301	Springhill Institution	NP Office	Welch Allyn	GS777	
69	Oxygen Concentrator		Springhill Institution		Philips	EverFlo Q	0135448
70	Oxygen Concentrator		Springhill Institution		Philips	Everflo Q	0056787
71	PHOROPTOR	BMD 11639	Springhill Institution	Optometrist	Topcon	238712	3617864
72	Portable Oxygen		Springhill Institution		OxyGo Portable Oxygen Concentrator	Inogen One G3 10-300	19492998
73	Printer-XRay	BMD 10414	Springhill Institution	X-ray	AGFA	Drystar 5300	
74	PROJECTOR-OPHTHALMIC	BMD 5348	Springhill Institution	Optometrist	Topcon	238723	262933G
75	PUMP-COMPRESSOR	BMD 7862	Springhill Institution	Treatment	Devilbiss	5650C Pulmoaide	C170892
76	PUMP-COMPRESSOR	BMD 7867	Springhill Institution	Treatment	Invacare Corp	IRC1196 - PASSPORT	59400000897
77	PUMP-COMPRESSOR	BMD 8081	Springhill Institution		Innovocorp	IRC1196	
78	PUMP-INFUSION		Springhill Institution		Abbott	Plum 360 Large volume	43457394
79	PUMP-INFUSION		Springhill Institution		Abbott	Plum 360 Large volume	43468679
80	PUMP-INFUSION		Springhill Institution		Abbott	Plum 360 Large volume	43468220
81	REFRIGERATOR-BIOLOGICAL	BMD 10351	Springhill Institution	Medication Room	Sanyo	MPR721	11110815
82	REFRIGERATOR-BIOLOGICAL	BMD 15198	Springhill Institution	Health Services	Panasonic	SANSR-L6111W	833-6-4219-527-000
83	REFRIGERATOR-BIOLOGICAL	BMD 8079	Springhill Institution	Treatment	Panasonic	SR-L6111W	120700510
84	REFRIGERATOR-BIOLOGICAL	BMD 8085	Springhill Institution	Treatment	Panasonic		121001043
85	SAW-CAST	BMD 10434	Springhill Institution	Treatment	BSN Medical Inc.	0295-200	28330
86	SCALE-PATIENT	BMD 10427	Springhill Institution		DET	6437	
87	SCALE-PATIENT	BMD 16032	Springhill Institution	hallway	Health O Meter	2500KL	2500L00415
88	STERILIZER-STEAM-TABLE-TOP	BMD 15199	Springhill Institution	Dentist	Scican	Statim 5000 G4	511017I00028

Annex M - Equipment List - Springhill Institution / Annexe M - Liste d'équipements - Établissement Springhill

	A	B	C	D	E	F	G
1	Item	Equipment Code Code d'équipement	Institution Établissement	Location Localisation	Manufacturer Nom du fabricant	Model Modèle	Serial Number Numéro de série
89	STRETCHER	BMD 10416	Springhill Institution	Treatment	Pedigo	PP7500-W-Spec	0199025
90	STRETCHER	BMD 16299	Springhill Institution	NP office	Stryker	CA LP 748	1612038569
91	STRETCHER	BMD 7861	Springhill Institution	Treatment	Hill Rom	268261	F093AW7612
92	STRETCHER		Springhill Institution	RRC	Stryker	CA LP 748	1612038570
93	TABLE-PHYSIO-TREATMENT	BMD 5335	Springhill Institution	Doctor's Office	Chattanooga Canada	TRT400	11717
94	TABLE-XRAY	BMD 11677	Springhill Institution	X-ray	GE	X171137 Picker	B25863
95	TONOMETER	BMD 16031	Springhill Institution	Optometrist	iCare	IC100	0598
96	Ultrasound-Diagnostic Bladder		Springhill Institution	Treatment	SeeMore	ViewBladder10	
97	Warmer-Blanket-Solution	BMD 7859	Springhill Institution	Treatment	Pedigo	P-2010S	545087
98	WHEELCHAIR	BMD 10503	Springhill Institution	Health Services	Cardinal Health	CW0006PEL	
99	WHEELCHAIR	BMD 10519	Springhill Institution	Health Services	AMG Sales Agency	Airgo Pro Care 700-631	
100	WHEELCHAIR	BMD 11958	Springhill Institution	Health Services	Cardinal Health	CW0071PEL	
101	WHEELCHAIR	BMD 11959	Springhill Institution	Health Services	Cardinal Health	CW0006PEL	
102	WHEELCHAIR	BMD 11960	Springhill Institution	Health Services	Cardinal Health	CW00071PEL	
103	WHEELCHAIR	BMD 16028	Springhill Institution	Health Services	Cardinal Health	CW0006PEL	
104	X-Ray Head	BMD 7849	Springhill Institution	dentist	Belmont	428384	EH0710201
105	XRAY-UNIT	BMD 11667	Springhill Institution	X-ray	GE	600-030	46335HL6
106	XRAY-UNIT	BMD 11681	Springhill Institution	X-ray	GE	5892919	4034
107	XRAY-UNIT	BMD 15902	Springhill Institution	X-ray	GENDEX	Orthoralix 8500	35-1858143DP
108	Zeiss Visuscout mobile Retinal Camera	TBD	Springhill Institution	TBD			